

Analytical Method Information

Analyte	MDL	Reporting Limit	Surrogate %R	Duplicate RPD	Matrix Spike %R	Matrix Spike RPD	Blank Spike / LCS %R	Blank Spike / LCS RPD
525.2 DHS SVOA Reg3 in Water (EPA 525.2)								
Preservation: Dechlor with 50mg/L Na2SO3, HCl to pH<2, <6°C								
Container: 1L Amber-HCl			Amount Required: 2 1L Ambers			Hold Time: 14 days		
Benzo (a) pyrene	0.10	0.10 ug/l			70 - 130	20	70 - 130	20
Di(2-ethylhexyl)adipate	5.0	5.0 ug/l			70 - 130	20	70 - 130	20
Di(2-ethylhexyl)phthalate	3.0	3.0 ug/l			70 - 130	20	70 - 130	20
surr: 1,3-Dimethyl-2-nitrobenzene				70 - 130				
surr: Perylene-d12				70 - 130				
surr: Triphenyl phosphate				70 - 130				